

# Cross-Image Component Linking Report

Image 1: ocr-rx1.png · Image 2: ocr-rx2.png

|                            | Image 1 | Image 2 |
|----------------------------|---------|---------|
| Components detected        | 4       | 3       |
| Equivalence classes        | 3       | 3       |
| Within-image relations     | 0       | 3       |
| Identity (very high) links | 0       |         |
| SKU (high) links           | 0       |         |
| Therapeutic (medium) links | 0       |         |
| Unmatched                  | 3       | 3       |

## Confidence tiers

|             |                                                                                                      |                  |
|-------------|------------------------------------------------------------------------------------------------------|------------------|
| Identity    | Same physical lot. Shared NDC + Lot (+ optional Exp), or shared Lot + Exp/Mfg with no conflicts.     | <i>very high</i> |
| SKU         | Same product SKU. Shared NDC; lots may differ or be unobserved.                                      | <i>high</i>      |
| Therapeutic | Same drug + strength but different manufacturer / NDC. Functionally substitutable, not the same SKU. | <i>medium</i>    |

Image 1 segmentation





Image 2 segmentation



## Within-image equivalence classes

Components grouped here share enough identifying fields (typically Lot + Exp, possibly with NDC or Mfg) to be treated as the same physical lot. Their evidence is merged before cross-image matching.

Image 1 · Components [2, 3] (grouped: same\_lot\_ndc)

|                                                                                   |              |                                          |                                                                                   |              |                                         |
|-----------------------------------------------------------------------------------|--------------|------------------------------------------|-----------------------------------------------------------------------------------|--------------|-----------------------------------------|
|  | Image / Obj  | Image 1 · Component 2                    |  | Image / Obj  | Image 1 · Component 3                   |
|                                                                                   | Type         | vial                                     |                                                                                   | Type         | vial                                    |
|                                                                                   | Manufacturer | —                                        |                                                                                   | Manufacturer | —                                       |
|                                                                                   | Product      | —                                        |                                                                                   | Product      | —                                       |
|                                                                                   | Strength     | —                                        |                                                                                   | Strength     | —                                       |
|                                                                                   | NDC          | —                                        |                                                                                   | NDC          | —                                       |
|                                                                                   | LOT          | K250250                                  |                                                                                   | LOT          | K250250                                 |
|                                                                                   | EXP          | 09/26                                    |                                                                                   | EXP          | 09/26                                   |
|                                                                                   | Volume       | 2 mL                                     |                                                                                   | Volume       | 2 mL                                    |
|                                                                                   | OCR          | rot 90° + -0.0° · PSM 3 · otsu · conf 92 |                                                                                   | OCR          | rot 0° + +2.0° · PSM 4 · otsu · conf 87 |

## Within-image relations

Components in the same image that aren't the same physical lot but represent the same SKU or are therapeutically equivalent.

### Image 2: [1] ↔ [2] sku (high)

shared: ndc, product, strength, volume\_ml



|              |                                         |
|--------------|-----------------------------------------|
| Image / Obj  | Image 2 · Component 1                   |
| Type         | bag                                     |
| Manufacturer | —                                       |
| Product      | Thiamine HCl Injection                  |
| Strength     | 200 mg/2 mL                             |
| NDC          | 04001-0243-01                           |
| LOT          | —                                       |
| EXP          | —                                       |
| Volume       | 2 mL                                    |
| OCR          | rot 0° + -1.0° · PSM 3 · otsu · conf 88 |

SKU



|              |                                          |
|--------------|------------------------------------------|
| Image / Obj  | Image 2 · Component 2                    |
| Type         | bag                                      |
| Manufacturer | —                                        |
| Product      | Thiamine HCl Injection                   |
| Strength     | 200 mg/2 mL                              |
| NDC          | 04001-0243-01                            |
| LOT          | —                                        |
| EXP          | —                                        |
| Volume       | 2 mL                                     |
| OCR          | rot 90° + -0.0° · PSM 3 · otsu · conf 91 |

### Image 2: [1] ↔ [3] therapeutic (medium)

shared: product, strength, volume\_ml

notes: same product/strength, different manufacturer



|              |                                         |
|--------------|-----------------------------------------|
| Image / Obj  | Image 2 · Component 1                   |
| Type         | bag                                     |
| Manufacturer | —                                       |
| Product      | Thiamine HCl Injection                  |
| Strength     | 200 mg/2 mL                             |
| NDC          | 04001-0243-01                           |
| LOT          | —                                       |
| EXP          | —                                       |
| Volume       | 2 mL                                    |
| OCR          | rot 0° + -1.0° · PSM 3 · otsu · conf 88 |

THERAPEUTIC



|              |                                          |
|--------------|------------------------------------------|
| Image / Obj  | Image 2 · Component 3                    |
| Type         | bag                                      |
| Manufacturer | Fresenius Kabi                           |
| Product      | Thiamine HCl Injection                   |
| Strength     | 200 mg/2 mL                              |
| NDC          | 03442-0500-10                            |
| LOT          | —                                        |
| EXP          | —                                        |
| Volume       | 2 mL                                     |
| OCR          | rot 0° + -1.0° · PSM 11 · otsu · conf 94 |

Image 2: [2] ↔ [3] therapeutic (medium)

shared: product, strength, volume\_ml

notes: same product/strength, different manufacturer



| Image / Obj  | Image 2 · Component 2                    |
|--------------|------------------------------------------|
| Type         | bag                                      |
| Manufacturer | —                                        |
| Product      | Thiamine HCl Injection                   |
| Strength     | 200 mg/2 mL                              |
| NDC          | 04001-0243-01                            |
| LOT          | —                                        |
| EXP          | —                                        |
| Volume       | 2 mL                                     |
| OCR          | rot 90° + -0.0° · PSM 3 · otsu · conf 91 |


THERAPEUTIC



| Image / Obj  | Image 2 · Component 3                    |
|--------------|------------------------------------------|
| Type         | bag                                      |
| Manufacturer | Fresenius Kabi                           |
| Product      | Thiamine HCl Injection                   |
| Strength     | 200 mg/2 mL                              |
| NDC          | 03442-0500-10                            |
| LOT          | —                                        |
| EXP          | —                                        |
| Volume       | 2 mL                                     |
| OCR          | rot 0° + -1.0° · PSM 11 · otsu · conf 94 |

Unmatched components


Image 1 · Component(s) [1] (no cross-image match — no\_candidate)



|              |                                              |
|--------------|----------------------------------------------|
| Image / Obj  | Image 1 · Component 1                        |
| Type         | bag                                          |
| Manufacturer | B. Braun                                     |
| Product      | 0.9% Sodium Chloride Injection, USP          |
| Strength     | 0.9 g                                        |
| NDC          | 00264-7800-10                                |
| LOT          | US021790                                     |
| EXP          | 04/27                                        |
| Volume       | 100 mL                                       |
| OCR          | rot 0° + -1.0° · PSM 11 · inv_otsu · conf 89 |


(none)  
ONLY  
IMG 1

Image 1 · Component(s) [2, 3] (no cross-image match — no\_candidate)



|              |                                          |
|--------------|------------------------------------------|
| Image / Obj  | Image 1 · Component 2                    |
| Type         | vial                                     |
| Manufacturer | —                                        |
| Product      | —                                        |
| Strength     | —                                        |
| NDC          | —                                        |
| LOT          | K250250                                  |
| EXP          | 09/26                                    |
| Volume       | 2 mL                                     |
| OCR          | rot 90° + -0.0° · PSM 3 · otsu · conf 92 |


(none)  
ONLY  
IMG 1



|              |                                         |
|--------------|-----------------------------------------|
| Image / Obj  | Image 1 · Component 3                   |
| Type         | vial                                    |
| Manufacturer | —                                       |
| Product      | —                                       |
| Strength     | —                                       |
| NDC          | —                                       |
| LOT          | K250250                                 |
| EXP          | 09/26                                   |
| Volume       | 2 mL                                    |
| OCR          | rot 0° + +2.0° · PSM 4 · otsu · conf 87 |



Image 1 · Component(s) [4] (no cross-image match — no\_candidate)



|              |                                          |
|--------------|------------------------------------------|
| Image / Obj  | Image 1 · Component 4                    |
| Type         | vial                                     |
| Manufacturer | —                                        |
| Product      | —                                        |
| Strength     | —                                        |
| NDC          | —                                        |
| LOT          | 306122                                   |
| EXP          | 08/31/2095                               |
| Volume       | 2 mL                                     |
| OCR          | rot 0° + -2.0° · PSM 11 · otsu · conf 82 |

ONLY  
IMG 1

(none)

Image 2 · Component(s) [1] (no cross-image match — no\_candidate)

(none)

ONLY  
IMG 2



|              |                                         |
|--------------|-----------------------------------------|
| Image / Obj  | Image 2 · Component 1                   |
| Type         | bag                                     |
| Manufacturer | —                                       |
| Product      | Thiamine HCl Injection                  |
| Strength     | 200 mg/2 mL                             |
| NDC          | 04001-0243-01                           |
| LOT          | —                                       |
| EXP          | —                                       |
| Volume       | 2 mL                                    |
| OCR          | rot 0° + -1.0° · PSM 3 · otsu · conf 88 |

Image 2 · Component(s) [2] (no cross-image match — no\_candidate)

(none)

ONLY  
IMG 2



|              |                                          |
|--------------|------------------------------------------|
| Image / Obj  | Image 2 · Component 2                    |
| Type         | bag                                      |
| Manufacturer | —                                        |
| Product      | Thiamine HCl Injection                   |
| Strength     | 200 mg/2 mL                              |
| NDC          | 04001-0243-01                            |
| LOT          | —                                        |
| EXP          | —                                        |
| Volume       | 2 mL                                     |
| OCR          | rot 90° + -0.0° · PSM 3 · otsu · conf 91 |

Image 2 · Component(s) [3] (no cross-image match — no\_candidate)

(none)

|               |                                                                                    |              |                                          |
|---------------|------------------------------------------------------------------------------------|--------------|------------------------------------------|
| ONLY<br>IMG 2 |  | Image / Obj  | Image 2 · Component 3                    |
|               |                                                                                    | Type         | bag                                      |
|               |                                                                                    | Manufacturer | Fresenius Kabi                           |
|               |                                                                                    | Product      | Thiamine HCl Injection                   |
|               |                                                                                    | Strength     | 200 mg/2 mL                              |
|               |                                                                                    | NDC          | 03442-0500-10                            |
|               |                                                                                    | LOT          | —                                        |
|               |                                                                                    | EXP          | —                                        |
|               |                                                                                    | Volume       | 2 mL                                     |
|               |                                                                                    | OCR          | rot 0° + -1.0° · PSM 11 · otsu · conf 94 |



## OCR appendix — raw text per component

### Image 1

**Component 1** (rot 0° + -1.0°, PSM 11, inv\_otsu, mean conf 89)

NDC 0264-7800-10 Baxter 0.9% Sodium Chloride Injection, USP 100 mL ~1 Each 100 mL contains Sodium Chloride 0.9g Water for Injection Qs. pH 5.0 (4.5 to 7.0) Calc. Osmolarity 308 mOsmoVL For intravenous use Single Dose Container Discard unused portion Store at 25°C (77°F); -3 excursions permitted to 15-30°C (59-86°F) [see USP Controlled Room Temperature]. LOT US021790 EXP 04/27 Baxter Healthcare Corporation Deerfield, IL 60015 USA Made in USA PL 146 Plastic 287123™ atk i

**Component 2** (rot 90° + -0.0°, PSM 3, otsu, mean conf 92)

LOT K250250 EXP 09/26 2mL

**Component 3** (rot 0° + +2.0°, PSM 4, otsu, mean conf 87)

LOT K250250 EXP 09/26 2mL

**Component 4** (rot 0° + -2.0°, PSM 11, otsu, mean conf 82)

LOT 306122% P 08/31/2025 2mL

### Image 2

**Component 1** (rot 0° + -1.0°, PSM 3, otsu, mean conf 88)

| Thiamine 200mg/2mL for Injection ua NDC 04001-0243-01

**Component 2** (rot 90° + -0.0°, PSM 3, otsu, mean conf 91)

Thiamine 200mg/2mL for Injection ad NDC 04001-0243-01

**Component 3** (rot 0° + -1.0°, PSM 11, otsu, mean conf 94)

Thiamine 200mg/2mL for Injection KABI FRESENIUS NDC 03442-0500-10